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**REF** Artikelnummer **LOT** Chargennummer

## Mouth Care Finger Stall – Instructions for Use

SILVER DUO SYSTEM: Textile finger stall with differently structured sides for cleaning and massaging remains antibacterial due to a silver finish.

Please read the 'safety information' and the 'information about application and care' thoroughly.

⚠ Please note instructions and packaging and keep it!

Made in Germany

#### 1. Scope of Application

These instructions for use apply to the following products:

**REF** 170-00 Mouth Care Finger Stall

#### 2. Contents

2 x mouth care finger stall, Instructions for use

#### 3. Purpose and intended Use

The mouth care finger stall is only intended for oral application. It's differently structured sides are suited for gum massages and cleaning the first minute teeth of babies. It can be applied to the first tooth and trains daily mouth care. Babies feel that gum massages are particularly soothing when the first minute teeth break through. Since the application can not be performed by the baby himself, the daily use strengthens the parent-child relationship. The mouth care finger stall can also be used by adults for cleaning the inside of the mouth and the tongue.

#### 4. Application and care instructions

Please clean the finger stall thoroughly with water approx. 40°C warm before the first use. Please put the finger stall over your index finger and moisten it. It is also suited for small amounts of tooth paste. The flexible material fits your finger perfectly without slipping. You can now control the pressure for a soothing massage (smooth side) or extensive cleaning (soft side) yourself and reach every part of the oral cavity. Rinse the finger stall under running warm water after every application. Then squeeze it well and leave to dry in a light, hygienic area which is open to the air. The practical double pack makes alternating use in the morning and evening possible. If you use it daily, replace it after 4 to 6 weeks. Please do not clean it in the washing machine or the dishwasher.

#### 5. Important safety information!

⚠ **For the safety of your child. CAUTION!**

No toy. Do not leave your child unattended with the finger stall. Check the finger stall before every single use. Replace it immediately, if you find the first signs of damage or defects!

#### 6. Material Composition

Fabric: 55 % polyester, 30 % polyamide, 15 % elastane

#### 7. Explanation of Symbols used

 Manufacturer  Attention!  Read Instructions for Use

**REF** Article Number **LOT** Batch number

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