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Ⓛ Artikelnummer Ⓛ Chargennummer

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### Ⓛ Natural finger toothbrush made from natural rubber – Instructions for Use

#### Dear Customer,

Thank you for buying this product. Below we have compiled important product information for you. Please read this carefully, in particular the safety advice, as well as the advice on how to use and care for the product.

⚠ Please follow the instructions and retain the instruction sheet and the packaging!

#### 1. Scope of Application

These instructions for use are applicable for the following

Ⓛ 687-00 Natural finger toothbrush made from natural rubber

#### 2. Contents

- Natural finger toothbrush made from natural rubber

- Instructions for Use

#### 3. Intended Purpose and Use

The finger toothbrush is designed to gently clean the teeth and massage the gums of babies and small children, and for supplementary teeth cleaning for children, adults and those in need of care. It can also be used to clean the tongue.

#### 4. Application

Pull the finger toothbrush over a suitable finger. Use the bristle side for cleaning the teeth and the dimpled, ridged side for massaging the gums and cleaning the tongue.

### IS THE FINGER TOOTHBRUSH TOO TIGHT?

Tips are available on our website!

#### 5. Important safety information!

##### ⚠ For the safety of your child CAUTION!

No toy. It must only be used by adults. Nursing staff should use it for those in need of care. Before using in the mouth, make sure that the finger toothbrush is firmly on your finger. Note that occasionally, an individual may have an allergic reaction to the natural rubber. Should you notice an allergic reaction while using the finger toothbrush, stop using it immediately. Check the finger toothbrush before and after use for damage. Dispose of the finger toothbrush at the first sign of damage, e.g. if there are cracks, changes in shape or the material is flaking off. Do not put the finger toothbrush in the freezer! Do not subject the finger toothbrush to direct sunlight, heat or continuous moisture. Avoid contact with greasy substances. Please note that when using toothpaste, the finger toothbrush may become slightly sticky. This is not a sign of an accelerated aging process, but is due to dental hygiene products (see the care and cleaning instructions).

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#### 6. Care and cleaning instructions

Clean the finger toothbrush thoroughly with warm running water before and after each use. Before using it for the first time, or if it is extremely dirty, it can be immersed in boiling water for approx. 30 seconds. When using toothpaste, make sure that you take extra care to remove all residue. Allow the finger toothbrush to air dry and store it in a dry, hygienic place. Do not use any abrasive cleaning agents or washing up liquid to clean the finger toothbrush. Not suitable for being heated up in a microwave or being cleaned in a dishwasher. The finger toothbrush may discolour from boiling or over time - this does not affect the quality of the product and is a natural process. Note that the finger toothbrush is made of natural rubber. The natural rubber smell will disappear with time and is completely harmless.

#### 7. Information about Shelf Life and Disposal

Replace the finger toothbrush every 4 weeks to ensure safety and hygiene. The finger toothbrush can be disposed of in household waste and recycled thermally without any problems.

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